


HU 6

Hydrolysis Unit

Convenient, reliable and safe hydrolysis unit,
as preliminary step for total fat determination.



Driven by curiosity



Acid & Base Hydrolysis of Fats

The Velp Hydrolysis Unit HU 6 is the ideal solution for sample preparation prior to solvent extraction for total fat determination according to official methods.

- Fast and precise heating with 0.5°C temperature stability and homogeneity
- Convenient operations saving operator time
- Faster and enhanced filtration thanks to the valves that allow to selectively interrupt the vacuum in each position
- Continuous filtration and rinsing of six samples in parallel
- Intuitive digital interface
- Low energy consumption
- Space saving design



Instrument - Codes

HU 6	230 V / 50-60 Hz	F30300150
HU 6	115 V / 50-60 Hz	F30310150

Technical Data

Construction material	Epoxy painted stainless steel structure
Number of samples	6 samples
Set temperature	Display visualization
Countdown	Display visualization
Display	LCD
Number of programs	20
Selectable language	IT, EN
Temperature range	From room temp. to 200 °C
Temperature calibration	Automatic
Stability of the heating block temperature	± 0,5°c
Precision of the heating block temperature	± 0,5°c
Reaching of set temperature	Acoustic and visual signal
Power	1150 W
Weight	14 kg / 30.8 lb
Dimensions (WxHxD)	420x230x130 mm (16.5x9.0x5.1 in)



Headquarters
Velp Scientifica Srl
Via Stazione 16
20865 Usmate (MB)
Italia
T +39 039 628811
velpitalia@velp.com

LATAM
velplatam@velp.com

USA
Velp Scientifica Inc
40, Burt Drive, Unit #1,
Deer Park
NY 11729 - U.S.
T +1 631 573 6002
velpusa@velp.com

SEA & Pacific
velpsea@velp.com

China
Velp China Co. Ltd.
Room 828, Building 1,
778 Jinji Road, Pudong,
Shanghai, China
Tel. +8621 34500630
velpchina@velp.com

India
Velp India Pvt. Ltd.
Office 323-324, Building 2
Raheja Tesla, MIDC, Turbhe,
400705, Navi Mumbai, India
velpindia@velp.com

Velp Official Partner

Rev 3.02.2026

